

# Product Characterization Sheet

## HH1067

Human Hepatocytes, Catalog Number 82001



### Classification

Plateability	Suspension
Number of days plateable	NA
Confluency	NA
P450 Inducibility	NA
Transporter activity	NA
Number of donors	1

### Donor Demographics

Gender	Male
Age	49
Race	Caucasian
Cause of death	Head Trauma
BMI	26.3
Smoking	No
Alcohol	Yes
Substance abuse	Yes
Medical history	None
Infectious diseases	HBV-, HCV-, HIV-, CMV-, EBV (IgG)+

### Post-thaw Viability and Yield

Viability	85 %
Yield	4.1 million

**Characterization:** Hepatocytes were thawed using 37°C UCRM™ and centrifuged for 10 minutes at 100g. After removing the supernatant, hepatocytes were re-suspended in 4°C HQM™ and counted for viability and yield using the Trypan Blue exclusion method.

### Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (μM)	Incubation Time (minutes)	Metabolite Quantified	Activity (pmol/minute/million cells)
CYP1A2	Phenacetin (100)	15	Acetaminophen	70.4 ± 11.4
CYP2A6	Coumarin (50)	30	7-Hydroxycoumarin	19.3 ± 0.8
CYP2B6	Bupropion (500)	15	Hydroxybupropion	7.9 ± 3.6
CYP2C8	Paclitaxel (20)	15	6α-Hydroxypaclitaxel	53.5 ± 20.6
CYP2C9	Diclofenac (25)	15	4-Hydroxydiclofenac	126 ± 5.4
CYP2C19	S-Mephenytoin (250)	30	4-Hydroxymephenytoin	6.7 ± 0.9
CYP2D6	Dextromethorphan (15)	15	Dextrorphan	14.5 ± 4.3
CYP2E1	Chlorzoxazone (250)	15	6-Hydroxychlorzoxazone	47.5 ± 4.9
CYP3A4	Midazolam (20)	10	1-Hydroxymidazolam	28 ± 4.3
	Testosterone (200)	15	6β-Hydroxytestosterone	461 ± 231
ECOD	7-Ethoxycoumarin (100)	30	7-Hydroxycoumarin	42.4 ± 8.9
UGT	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin glucuronide	474 ± 60.4
Sulfotransferase	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin sulfate	11.5 ± 0.9

**CYP450 Activity Assessment:** The hepatocytes were incubated at a cell density of 0.5 million cells/mL in a 48-well plate (125,000 hepatocytes/well) for the designated time durations with isoform-selective substrates. The metabolites were identified and analyzed using LC-MS/MS.

IVAL cell culture media and tissue culture plates used in this evaluation:

- Recovery of thawed hepatocytes - Cat. No. 81015 - UCRM™ Universal Cryopreservation Recovery Media, 50 mL tube
- Suspension and incubation of hepatocytes - Cat. No. 81039/81040 - HQM™ Hepatocyte Incubation Media, 50 mL tube/500 mL bottle

To inquire about our products and services or for technical questions please contact:

- In Vitro ADMET Laboratories by phone at +1 (866) 458-1094 or +1 (410) 869-9037 or email at [info@invitroadmet.com](mailto:info@invitroadmet.com)