

# Product Characterization Sheet

## HH1105

Human Hepatocytes, Catalog Number 82002



### Classification

Plateability	Suspension
Number of days plateable	NA
Confluency	NA
P450 Inducibility	NA
Transporter activity	NA
Number of donors	1

### Donor Demographics

Gender	Female
Age	59 years
Race	Asian
Cause of death	CVA
BMI	21.4
Smoking	No
Alcohol	No
Substance abuse	No
Medical history	Diabetes, HTN
Infectious diseases	HBV-, HCV-, HIV-, CMV+, EBV (IgG)+

### Post-thaw Viability and Yield

Viability	83 %
Yield	5.6 million

**Characterization:** Hepatocytes were thawed using 37°C UCRM™ and centrifuged for 10 minutes at 100g. After removing the supernatant, hepatocytes were re-suspended in 4°C HQM™ and counted for viability and yield using the Trypan Blue exclusion method.

### Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (μM)	Incubation Time (minutes)	Metabolite Quantified	Activity (pmol/minute/million cells)
CYP1A2	Phenacetin (100)	15	Acetaminophen	10.4 ± 1.3
CYP2A6	Coumarin (50)	30	7-Hydroxycoumarin	26 ± 1.6
CYP2B6	Bupropion (500)	15	Hydroxybupropion	0.84 ± 0.73
CYP2C8	Paclitaxel (20)	15	6α-Hydroxypaclitaxel	0 ± 0.2
CYP2C9	Diclofenac (25)	15	4-Hydroxydiclofenac	63.9 ± 14.4
CYP2C19	S-Mephenytoin (250)	30	4-Hydroxymephenytoin	0 ± 0
CYP2D6	Dextromethorphan (15)	15	Dextrorphan	3.13 ± 1.24
CYP2E1	Chlorzoxazone (250)	15	6-Hydroxychlorzoxazone	15.7 ± 2.5
CYP3A4	Midazolam (20)	10	1-Hydroxymidazolam	42.3 ± 2.8
	Testosterone (200)	15	6β-Hydroxytestosterone	165 ± 17.6
ECOD	7-Ethoxycoumarin (100)	30	7-Hydroxycoumarin	10.6 ± 1.7
UGT	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin glucuronide	543 ± 33
Sulfotransferase	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin sulfate	10.8 ± 1.7

**CYP450 Activity Assessment:** The hepatocytes were incubated at a cell density of 0.5 million cells/mL in a 48-well plate (125,000 hepatocytes/well) for the designated time durations with isoform-selective substrates. The metabolites were identified and analyzed using LC-MS/MS.

IVAL cell culture media and tissue culture plates used in this evaluation:

- Recovery of thawed hepatocytes - Cat. No. 81015 - UCRM™ Universal Cryopreservation Recovery Media, 50 mL tube
- Suspension and incubation of hepatocytes - Cat. No. 81039/81040 - HQM™ Hepatocyte Incubation Media, 50 mL tube/500 mL bottle

To inquire about our products and services or for technical questions please contact:

- In Vitro ADMET Laboratories by phone at +1 (866) 458-1094 or +1 (410) 869-9037 or email at [info@invitroadmet.com](mailto:info@invitroadmet.com)