

Product Characterization Sheet

HH1037

Human Hepatocytes, Catalog Number 82001



Classification

Plateability	Suspension
Number of days plateable	NA
Confluency	NA
P450 Inducibility	NA
Transporter activity	NA
Number of donors	1

Donor Demographics

Gender	Male
Age	45 years
Race	Caucasian
Cause of death	CVA 2 nd to ICH
BMI	25.5
Smoking	Yes
Alcohol	Yes
Substance abuse	Yes
Medical history	HTN
Infectious diseases	HBV-, HCV-, HIV-, CMV+

Post-thaw Viability and Yield

Viability	88 %
Yield	3.3 million

Characterization: Hepatocytes were thawed using 37°C UCRM™ and centrifuged for 10 minutes at 100g. After removing the supernatant, hepatocytes were re-suspended in 4°C HQM™ and counted for viability and yield using the Trypan Blue exclusion method.

Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (μM)	Incubation Time (minutes)	Metabolite Quantified	Activity (pmol/minute/million cells)
CYP1A2	Phenacetin (100)	15	Acetaminophen	83.7
CYP2A6	Coumarin (50)	30	7-Hydroxycoumarin	45.2
CYP2B6	Bupropion (500)	15	Hydroxybupropion	6.8
CYP2C8	Paclitaxel (20)	15	6α-Hydroxypaclitaxel	8.2
CYP2C9	Diclofenac (25)	15	4-Hydroxydiclofenac	20.6
CYP2C19	S-Mephenytoin (250)	30	4-Hydroxymephenytoin	9.7
CYP2D6	Dextromethorphan (15)	15	Dextrorphan	5.2
CYP2E1	Chlorzoxazone (250)	15	6-Hydroxychlorzoxazone	9.4
CYP3A4	Midazolam (20)	15	1-Hydroxymidazolam	23.4
	Testosterone (200)	10	6β-Hydroxytestosterone	256.7
ECOD	7-Ethoxycoumarin (100)	30	7-Hydroxycoumarin	67.8
UGT	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin glucuronide	74.0
Sulfotransferase	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin sulfate	24.5

CYP450 Activity Assessment: The hepatocytes were incubated at a cell density of 0.5 million cells/mL in a 48-well plate (125,000 hepatocytes/well) for the designated time durations with isoform-selective substrates. The metabolites were identified and analyzed using LC-MS/MS.

IVAL cell culture media and tissue culture plates used in this evaluation:

- Recovery of thawed hepatocytes - Cat. No. 81015 - UCRM™ Universal Cryopreservation Recovery Media, 50 mL tube
- Suspension and incubation of hepatocytes - Cat. No. 81039/81040 - HQM™ Hepatocyte Incubation Media, 50 mL tube/500 mL bottle

To inquire about our products and services or for technical questions please contact:

- In Vitro ADMET Laboratories by phone at +1 (866) 458-1094 or +1 (410) 869-9037 or email at info@invitroadmet.com