

# Product Characterization Sheet

## 032-237

Rat Hepatocytes, Catalog Number 82019



### Classification

Plateability	Plateable
Number of days plateable	5
Confluency	100 %
P450 Inducibility	NA
Number of donors	4

### Lot Demographics

Gender	Female
Strain	Sprague Dawley
Age	7-10 weeks
Weight	200-250 g
Infectious diseases	Negative

### Post-thaw Viability and Yield

Viability	93 %
Yield	5.8 million

**Characterization:** Hepatocytes were thawed using 37°C UCRM™ and centrifuged for 5 minutes at 100g. After removing the supernatant, hepatocytes were re-suspended in UPCM™ and counted for viability and yield using the Trypan Blue exclusion method. Cells were plated at a density of 0.7 million cells per mL in a hand-coated collagen plate and allowed to attach 4-6 hours prior to a Matrigel® overlay.


### Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (µM)	Incubation Time (minutes)	Metabolite Quantified	Activity (pmol/minute/million cells)
ECOD	7-Ethoxycoumarin (100)	30	7-Hydroxycoumarin	154.70 ± 6.80
UGT	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin glucuronide	1288.54 ± 59.96
Sulfotransferase	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin sulfate	67.37 ± 3.99


**CYP450 Activity Assessment:** The hepatocytes were incubated at a cell density of 0.5 million hepatocytes/mL in a 12-well plate (500,000 hepatocytes/well) for the designated time durations with isoform-selective substrates. The metabolites were identified and analyzed using LC-MS/MS.

### Photomicrographs (100X, Phase Contrast)

Phase Contrast Day 2



Phase Contrast Day 4



**Monolayer Comments:** 032-237 has a high attachment efficiency at 90 % and continues to develop a monolayer confluency of 100 % by day 3. This lot exhibits excellent morphology and remains intact for more than 5 days in culture.

IVAL cell culture media and tissue culture plates used in this evaluation:

- Recovery of thawed hepatocytes - Cat. No. 81015 - UCRM™ Universal Cryopreservation Recovery Media, 50 mL tube
- Initial plating of hepatocytes - Cat. No. 81016 – UPCM™ IVAL Universal Primary Cell Plating Medium, 50 mL tube
- Sandwich culture with 0.25 mg Matrigel® - Cat. No. 81037/81038 - RHIM™ Rodent Hepatocyte Induction Media, 50 mL tube/500 mL bottle
- Suspension and incubation of hepatocytes - Cat. No. 81039/81040 - HQM™ Hepatocyte Incubation Media, 50 mL tube/500 mL bottle
- Collagen coated plates - Cat. No. 71006 - CellAffix™ 24-well Collagen Hand Coated tissue culture plate, 5 plates per pack

To inquire about our products and services or for technical questions please contact:

- In Vitro ADMET Laboratories by phone at +1 (866) 458-1094 or +1 (410) 869-9037 or email at [info@invitroadmet.com](mailto:info@invitroadmet.com)