

Product Characterization Sheet

PCEX3096

MetMax™ Monkey Enterocytes, Catalog Number 82141



Specification

Number of donors	3
Volume	2.5 mL
Cell concentration	2.0 x 10 ⁶ /mL
Plateability	Suspension

Demographics

Gender	Female
Strain	Cynomolgus
Age	3-4 years
Infectious diseases	Negative

Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (μM)	Metabolite Quantified	Activity (pmol/minute/ million cells)
CYP1A2	Phenacetin (100)	Acetaminophen	86.2 ± 9.34
CYP2A6	Coumarin (50)	7-HC, 7-HC-Sulfate, 7-HC-Glucuronide	0 ± 0
CYP2A6	Coumarin (50)	7-Hydroxycoumarin	0 ± 0
CYP2A6	Coumarin (50)	7-Hydroxycoumarin Glucuronide	0 ± 0
CYP2A6	Coumarin (50)	7-Hydroxycoumarin Sulfate	0 ± 0
CYP2B6	Bupropion (500)	Hydroxybupropion	0.129 ± 0.2
CYP2C8	Paclitaxel (20)	6α-hydroxypaclitaxel	0.02 ± 0.01
CYP2C9	Diclofenac (25)	4-OH Diclofenac	2.53 ± 0.1
CYP2C19	S-Mephenytoin (250)	4-OH S-Mephenytoin	0 ± 0
CYP2D6	Dextromethorphan (15)	Dextrophan	0 ± 0
CYP2E1	Chlorzoxazone (250)	6-OH Chlorzoxazone	3.59 ± 0.05
CYP3A4-1	Midazolam (20)	1-Hydroxymidazolam	1.81 ± 0.3
CYP3A4-2	Testosterone (200)	6β-hydroxytestosterone	0 ± 0
ECOD	7-Ethoxycoumarin (100)	7-HC, 7-HC-Sulfate, 7-HC-Glucuronide	398 ± 13.9
ECOD	7-Ethoxycoumarin (100)	7-Hydroxycoumarin	398 ± 13.5
ECOD	7-Ethoxycoumarin (100)	7-Hydroxycoumarin Glucuronide	0.66 ± 0.12
ECOD	7-Ethoxycoumarin (100)	7-Hydroxycoumarin Sulfate	0.4 ± 0.09
UGT	7-Hydroxycoumarin (100)	7-Hydroxycoumarin Glucuronide	4.72 ± 0.43
SULT	7-Hydroxycoumarin (100)	7-Hydroxycoumarin Sulfate	1.5 ± 0.23
GST	Acetaminophen (10)	Acetaminophen Glutathione	1.4 ± 0.8
UGT	Acetaminophen (10)	Acetaminophen Glucuronide	0.13 ± 0.1
SULT	Acetaminophen (10)	Acetaminophen Sulfate	0 ± 0
FMO	Benzydamine HCl (250)	Benzydamine-N-Oxide	2.51 ± 0.2
MAO	Kynuramine HCl (160)	4-hydroxyquinoline	14.2 ± 4.9
AO	Carbazeran HCl (10)	4-Hydroxycarbazeran	0.033 ± 0.01
NAT1	4-Aminobenzoic HCl (200)	N-Acetyl-p-aminobenzoic acid	17.6 ± 1.12
NAT2	Sulfamethazine (100)	N-Acetyl-sulfamethazine	2.18 ± 0.17
1A1	7-EROD (100)	Resozufin	0 ± 0
2J2	Astemizole (50)	O-Demethyl Astemizole	4.31 ± 0.13
CES2	Irinotecan (50)	SN38	1.05 ± 0.08

CYP450 Activity Assessment: MetMax™ enterocytes were thawed in a 37°C water bath. A 1000 μL Eppendorf pipette (with tip) was used to pipet up and down three times to achieve homogeneity. MetMax™ enterocytes were prepared at a density of 2 million cells per mL. Equal volumes of MetMax™ enterocytes and 2X substrate were incubated for 15 minutes with isoform-selective substrates at 37°C. The final cell density of MetMax™ enterocytes and substrate was 1.0 million cells/mL in a 96-well plate (96 well plates pre-loaded with 50 μL of 2X substrates per well, and metabolism initiated by adding 50 μL of MetMax™ enterocytes to each of the wells, 100,000 hepatocytes/well). The metabolites were identified and analyzed using LC-MS/MS.

IVAL cell culture media and tissue culture plates used in this evaluation:

- Substrate preparation for incubation of hepatocytes - Cat. No. 81039/81040 - HQM™ Hepatocyte Incubation Media, 50 mL tube/500 mL bottle

To inquire about our products and services or for technical questions please contact:

- In Vitro ADMET Laboratories by phone at +1 (866) 458-1094 or +1 (410) 869-9037 or email at info@invitroadmet.com

In Vitro ADMET Laboratories | Columbia, MD 410-869-9037 | Malden, MA 781-397-9300 | info@invitroadmet.com | www.invitroadmet.com