

Product Characterization Sheet

HH1029

Human Hepatocytes, Catalog Number 82001



Classification

Plateability	Suspension
Number of days plateable	NA
Confluency	NA
P450 Inducibility	NA
Transporter activity	NA
Number of donors	1

Donor Demographics

Gender	Male
Age	51 years
Race	Caucasian
Cause of death	Anoxia 2 nd to CVA
BMI	18.0
Smoking	Yes
Alcohol	Yes
Substance abuse	Yes
Medical history	Arrhythmia
Infectious diseases	HBV-, HCV-, HIV-, CMV+

Post-thaw Viability and Yield

Viability	84 %
Yield	4.5 million

Characterization: Hepatocytes were thawed using 37°C UCRM™ and centrifuged for 10 minutes at 100g. After removing the supernatant, hepatocytes were re-suspended in 4°C HQM™ and counted for viability and yield using the Trypan Blue exclusion method.

Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (µM)	Incubation Time (minutes)	Metabolite Quantified	Activity (pmol/minute/million cells)
CYP1A2	Phenacetin (100)	15	Acetaminophen	57.0
CYP2A6	Coumarin (50)	30	7-Hydroxycoumarin	110.9
CYP2B6	Bupropion (500)	15	Hydroxybupropion	24.6
CYP2C8	Paclitaxel (20)	15	6α-Hydroxypaclitaxel	0.2
CYP2C9	Diclofenac (25)	15	4-Hydroxydiclofenac	51.0
CYP2C19	S-Mephenytoin (250)	30	4-Hydroxymephenytoin	1.0
CYP2D6	Dextromethorphan (15)	15	Dextrorphan	90.8
CYP2E1	Chlorzoxazone (250)	15	6-Hydroxychlorzoxazone	285.7
CYP3A4	Midazolam (20)	10	1-Hydroxymidazolam	8.8
	Testosterone (200)	15	6β-Hydroxytestosterone	81.9
ECOD	7-Ethoxycoumarin (100)	30	7-Hydroxycoumarin	4.9
UGT	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin glucuronide	397.7
Sulfotransferase	7-Hydroxycoumarin (100)	30	7-Hydroxycoumarin sulfate	10.3

CYP450 Activity Assessment: The hepatocytes were incubated at a cell density of 0.5 million cells/mL in a 48-well plate (125,000 hepatocytes/well) for the designated time durations with isoform-selective substrates. The metabolites were identified and analyzed using LC-MS/MS.

IVAL cell culture media and tissue culture plates used in this evaluation:

- Recovery of thawed hepatocytes - Cat. No. 81015 - UCRM™ Universal Cryopreservation Recovery Media, 50 mL tube
- Suspension and incubation of hepatocytes - Cat. No. 81039/81040 - HQM™ Hepatocyte Incubation Media, 50 mL tube/500 mL bottle

To inquire about our products and services or for technical questions please contact:

- In Vitro ADMET Laboratories by phone at +1 (866) 458-1094 or +1 (410) 869-9037 or email at info@invitroadmet.com