

# Product Characterization Sheet

## HE3065

Human Enterocytes, Catalog Number 82097



### Classification

Plateability	Suspension
Number of days plateable	NA
Confluency	NA
P450 Inducibility	NA
Number of donors	1
Enterocyte size	13.9 ± 2.1

### Donor Demographics

Gender	Female
Age	Gestational period of 39.3 weeks
Race	Caucasian
Cause of death	Anencephaly
BMI	NA
Smoking	No
Alcohol	No
Substance abuse	No
Medical history	NA
Infectious diseases	NA

### Post-thaw Viability and Yield

Viability	79 %
Yield	0.75 million

**Characterization:** Enterocytes were thawed using 37°C CERM™ and centrifuged for 10 minutes at 100g. After removing the supernatant, enterocytes were re-suspended in 4°C HQM™ and counted for viability and yield using the Trypan Blue exclusion method.


HE3065 is the same donor as hepatocyte lot HH1060.

### Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (μM)	Incubation Time (minutes)	Metabolite Quantified	Activity (pmol/minute/million cells)
CYP2C9	Diclofenac (25)	60	4-Hydroxydiclofenac	0 ± 0
CYP2C19	S-Mephenytoin (250)	60	4'-Hydroxymephenytoin	0.01 ± 0.02
CYP3A4	Testosterone (200)	60	6β-Hydroxytestosterone	2.08 ± 0.32
	Midazolam (20)	60	1-Hydroxymidazolam	0.01 ± 0
UGT	7-Hydroxycoumarin (100)	60	7-Hydroxycoumarin glucuronide	-0.1 ± 0
Sulfotransferase	7-Hydroxycoumarin (100)	60	7-Hydroxycoumarin sulfate	-0.02 ± 0
CYP2J2	Astemizole (50)	60	O-Desmethyl astemizole	0.07 ± 0.13
CES2	Irinotecan (50)	60	SN38	0.19 ± 0.23

**CYP450 Activity Assessment:** The enterocytes were incubated at a cell density of 1.5 million enterocytes/mL in a 96-well plate (150,000 enterocytes/well) for the designated time duration of 120 minutes with isoform-selective substrates. The metabolites were identified and analyzed using LC-MS/MS.

### Photomicrographs (100X, Brightfield)



**Enterocyte Morphology:** Enterocyte cell size is approximately 12 – 17 microns in diameter. Lymphocytes are approximately 4 microns and may also be observed in the enterocyte cell population.

IVAL cell culture media and tissue culture plates used in this evaluation:

- Recovery of thawed enterocytes - Cat. No. 81056 - CERM™ Cryopreserved Enterocyte Recovery Media, 50 mL tube
- Suspension and incubation of enterocytes - Cat. No. 81039 - HQM™ Hepatocyte and Enterocyte Incubation Media, 50 mL tube

To inquire about our products and services or for technical questions please contact:

- In Vitro ADMET Laboratories by phone at +1 (866) 458-1094 or +1 (410) 869-9037 or email at [info@invitroadmet.com](mailto:info@invitroadmet.com)